

# ioGlutamatergic Neurons

Digital flyer  
Cat.No io1001  
(formerly e001)

# Discovery accelerated with reliable, scalable and functional human excitatory neurons


ioGlutamatergic Neurons have been reprogrammed from human induced pluripotent stem cells (iPSCs) using our precise reprogramming technology: opti-ox™<sup>1</sup> (optimized inducible overexpression). Human stem cells, within days, convert into consistent, mature, functional glutamatergic neurons providing a high quality human model for the study of neurological activity and disease.

1. Pawlowski et al. Stem Cell Reports 2017  
[www.ncbi.nlm.nih.gov/pmc/articles/PMC5390118](http://www.ncbi.nlm.nih.gov/pmc/articles/PMC5390118)

## Characteristics that inspire confidence



**Defined**  
Characterised by ICC and gene expression.



**Scalable**  
Industrial scale quantities at a price point that allows the cells to be used from research to screening scale.



**Consistent**  
Batch to batch reproducibility and homogeneity create a stable human model for excitatory neuronal activity and disease.



**Quick**  
Spontaneous electrical activity observed within 8 days post revival, and functional neuronal networks observed within 13 days.




**Easy-to-use**  
Cells arrive programmed to rapidly mature upon revival. One medium required in a two-step protocol.



**Cost-effective**  
Available in two vial sizes, tailored to suit your experimental needs with minimal waste.







Neuronal networks within two weeks post-revival.  
Immunocytochemistry demonstrating expression  
of TUBB3 at Day 14 post-revival.

10x; scale bar: 500µm  
Courtesy Paul Holloway, University of Oxford





# Mature neuronal properties within two weeks of induction

**Express glutamatergic neuron-specific markers & exhibit neurite outgrowth**

ioGlutamatergic Neurons are derived from human iPSCs by Neurogenin-2 (NEUROG2) driven opti-ox™ reprogramming, resulting in the rapid differentiation of stem cells into mature and functional neurons.

Cultures consist mainly of glutamatergic neurons (>80%) characterised by homogenous expression of pan-neuronal proteins MAP2 and TUBB3, and glutamatergic-specific transporter VGLUT2. The minor remaining fraction of the neuronal population express marker genes of cholinergic neurons.


A bulk RNA-seq analysis shows that ioGlutamatergic Neurons have a rostral CNS identity and express the classical cortical marker genes FOXP1 and TBR1 (data not shown).

**“ These cells provide a reliable and pure source of glutamatergic neurons, resembling primary human ones. They are ready-to-use which makes it much more easy for tissue culture work and for reproducible results.”**

Dr Koby Baranes, University of Cambridge



Immunocytochemistry at Day 11 post-revival





# Robust and scalable cells for high-throughput applications

**Validated for CellTiter-Glo® and HTRF® assays**

Cytotoxicity CellTiter-Glo® (CTG) and TR-FRET (HTRF®) assays for AKT serine/threonine kinase 1 (AKT) and Huntingtin (HTT) proteins were performed on ioGlutamatergic Neurons in 384-well plates treated with tool compound (cmp) at day 9 post-revival. Compound titration results in a concentration response curve for all three assays (mean±sd of 2 replicates).

CTG assay on ioGlutamatergic Neurons shows an excellent average signal/background ratio and high suitability for HTS. HTRF® assays on ioGlutamatergic Neurons show lower signals but with low variability, and could therefore also provide a suitable platform for HTS.

Iovino M., et al., 2019, Charles River Laboratories



Signal/Background		
CTG	CellTiter-Glo®	327
HTT	HTT HTRF®	1.65
AKT	AKT HTRF®	1.95



# Neuronal activity that matures over-time


## Examples of MaxOne high-resolution multi-electrode array (MEA) recordings of ioGlutamatergic Neurons in BrainPhys™ media

The activity maps show firing rate (A), spike amplitude (B) and % of active electrodes (C). Results demonstrate a time-dependent increase of spontaneous activity during neuronal maturation from 2 to 3 weeks post-revival.

Iovino M., et al., 2019, Charles River Laboratories

“These cells enable us to move rapidly as from the moment of plating within 4–7 days we have mature and functional neurons.”

Dr Shushant Jain, Charles River






# Rapid gain of functional activity

## Time-dependent increase of spontaneous bursting activity over a three week period

- A. The graph shows the % of active bursting electrodes for each time point.
- B. An example of a spontaneous spike, taken at Day 8 post-revival (1 second sweep, 32  $\mu$ V/-18  $\mu$ V).
- C. An example of a bursting phenotype, taken at Day 20 post-revival (1 second sweep, 16  $\mu$ V/-16  $\mu$ V). Cells were cultured in the standard user manual medium (Complete Glutamatergic Neuron Medium) and recorded on 64-electrode MEA (Multi Channel Systems).

NDimension (Science and Engineering) Ltd






ioGlutamatergic  
Neurons co-cultured  
with rat-derived  
astrocytes  
demonstrate  
a time-dependent  
increase in  
synchronous activity



Array Wide Spike Detection Rate histograms (AWSDR – a graphical measure of synchrony) for 10-minute recordings on Day 8, 13 and 20 post-revival, in co-culture with primary rat-derived astrocytes

Results show prominent synchronicity on Day 13, exemplified by the ‘spikier’ nature of the associated AWSDR, which increases at Day 20. Cells were cultured in the standard user manual medium (Complete Glutamatergic Neuron Medium) and recorded on 64-electrode MEA (Multi Channel Systems).

NDimension (Science and Engineering) Ltd





# Ready-to-culture and easy-to-use


ioGlutamatergic Neurons are delivered in a cryopreserved format and are programmed to rapidly mature upon revival

Cells are kept in a single culture medium formulation from revival with fully disclosed composition allowing modifications to fit customer’s bespoke experiments. The protocol for the generation of these cells is a three-phase process:

- 1. Induction (carried out at bit.bio)
- 2. Stabilization for 4 days with Doxycycline
- 3. Maintenance during which the neurons mature

“These cells were really easy-to-use and ready to be plated in the assay format.”

Dr Mariangela Iovino, Charles River



# Cost effective and flexible



ioGlutamatergic Neurons are compatible with plates ranging from 6 to 384 wells and are available in two vial sizes, tailored to suit your experimental needs with minimal waste

Recommended seeding density is 30,000 cells/cm<sup>2</sup>.

One Small vial can plate a minimum of 0.5 x 24-well plate, 0.75 x 96-well plate, or 1 x 384-well plate.



One Large vial can plate a minimum of 1 x 24-well plate, 1.5 x 96-well plate, or 2 x 384-well plates.

Recommended seeding density

30,000 cells/cm<sup>2</sup>



Small vial size

0.75 x 10<sup>6</sup> cells



Large vial size

1.5 x 10<sup>6</sup> cells





# Upcoming products

Our range is constantly expanding. Coming soon are iPSC-derived inhibitory ioGABAergic Neurons and neuronal cells with disease-relevant mutations for modelling Huntington’s, Parkinson’s, Frontotemporal Dementia and ALS.

Request an invite to be the first to know about our latest products by emailing [info@bit.bio](mailto:info@bit.bio).

# Contact us

To order or speak with a member of our team, email [orders@bit.bio](mailto:orders@bit.bio).

bit.bio  
The Dorothy Hodgkin Building  
Babraham Research Campus  
Cambridge CB22 3FH  
United Kingdom  
+44 (0) 1223 787 297

[www.bit.bio](http://www.bit.bio)



# Product specifications

<b>Starting material</b> Human iPSC line	<b>Recommended seeding density</b> 30,000 cells/cm <sup>2</sup>
<b>Donor</b> Caucasian adult male (skin fibroblast)	<b>Seeding compatibility</b> 6 to 384 well plates
<b>Differentiation method</b> opti-ox™ cellular reprogramming	<b>Quality control</b> Sterility, protein expression (IF) and gene expression (RT-qPCR)
<b>Karyotype</b> Normal (46, XY)	<b>User storage</b> LN2 or -150°C
<b>Vial size</b> Small: >0.75 × 10 <sup>6</sup> viable cells Large: >1.5 × 10 <sup>6</sup> viable cells	<b>Shipping info</b> Dry ice
	<b>Product use</b> These cells are for research use only

# Product applications

<ul style="list-style-type: none"><li>• Academic research</li><li>• Drug development</li><li>• Neurotoxicology</li><li>• High-throughput screening</li><li>• Genetic screening (e.g. CRISPR screening)</li></ul>	<b>Validated techniques include</b> <ul style="list-style-type: none"><li>• HTRF® and CellTiter-Glo® assays</li><li>• RT-qPCR and branch DNA</li><li>• Immunocytochemistry assays (ICC)</li><li>• Patch clamp and multi-electrode array (MEA)</li><li>• Bulk &amp; single cell RNA seq</li><li>• Co-culture</li><li>• 3D bioprinting</li></ul> <p>HTRF® is the registered trademark of Cisbio. CellTiter-Glo® is the registered trademark of Promega.</p>
--	---

# Pricing

Size	Price	
Small	£1,199.00	
Large	£1,799.00	
Multipack (4 Large)	£5,600.00	Please contact us for academic and non-profit discount

