


Ronco High Security Solution for Responder 5 Nurse Call System

High Security Key Lock Station R5KLS-4



The Key Lock Station (R5KLS-4) is a flush-mounted unit for use in psychiatric wards and other high security hospital areas. It is generally positioned outside the entrance to a room or ward and is used to activate a dry contact pushbutton located in the room. The R5KLS has four separate inputs and outputs and can accommodate up to four separate call priorities from a room. All call placed via the R5KLS-4 must be cancelled at the R5KLS-4. A built-in “LED” indicator will illuminate when the key is to the “on” position, activating the R5KL. When a call is placed, a second built-in “LED” will illuminate at the key lock station. The call will be placed through the Responder 5 system, annunciating the call on the nurse call console and the call systems associated equipment. The call is terminated when the staff cancels the call from the R5KLS-4 station by turning the key lock to the “off” position.

The R5KLS-4 is built on a .050” stainless steel plate and incorporates circuitry for continuously reliable operation and requires a separate 15VDC power supply, not included. See typical installation drawings for details.

Dimensions: standard single gang plate: 2.75”W x 4.5”H x 1.5” D


Voltage Requirements: 12- 15.5 VDC, 50 ma all four outputs active

Current Draw @ 15VDC: 10 MA on/standby, 50 MA all four outputs active.

Back Box: fits standard single gang electrical box and single gang cut-in box.

R5KLS-4 input works with dry contact closures only and used with R5HSSCP-1 and R5WPPB.

Ronco High Security Solution for Responder 5 Nurse Call System




Ronco High Security Solution for Responder 5 Nurse Call System

High Security Staff Station Cover Plate R5HSSCP-1




The High Security Steel Cover Plate (R5HSSCP-1) is a 3 gang, .050", stainless steel plate, designed to protect the R5 Staff Station (353101) when used in psychiatric wards and other high security hospital areas. The plate incorporates a vandal resistant dry contact switch which is used in conjunction with the R5KLS and is used to initiate a call from the room. The plate installs over the R5 Staff Station (353101) without drilling or modification to the station.

Dimension: 6.54"W x 4.69"H




Ronco High Security Solution for Responder 5 Nurse Call System

High Security Push Button Cover Plate R5HSSCP-2



The High Security Steel Cover Plate (R5HSSCP-2) is a 3 gang, .050", stainless steel plate, designed to protect the R5 Staff Station when used in psychiatric wards and other high security hospital areas. The plate installs over the R5 Staff Station (353101) without modification to the station. The R5HSSCP-2 allows access and operation of the station's call and cancel switch. The R5KLS-4 is not used with this plate configuration.

Dimension: 6.54"W x 4.69"H




Ronco High Security Solution for Responder 5 Nurse Call System

High Security Push Button Cover Plate R5HSPBP



The High Security Push Button Cover Plate (R5HSPBCP) is a .050", stainless steel plate, designed to protect Rauland R5 Code (354011) and Staff Assist Station (354012) when used in psychiatric wards and other high security hospital areas. The plate installs over the push button stations without drilling or modification to the station. This plate is NOT USED with the R5KLS-4.



Dimensions: Cover Plate: 4.49"W x 4.70"H
Back Plate: 4.25"W x 4.25"H


Back Box: single gang masonry style (will not fit ganged cut-in box), 2 ½ deep

Ronco High Security Solution for Responder 5 Nurse Call System

High Security Cancel Station Cover Plate R5HSCSCP



The High Security Cancel Station Cover Plate (R5HSCSCP) is a .050", stainless steel plate, designed to protect Rauland R5 Cancel Station (354010) when used in psychiatric wards and other high security hospital areas. The plate installs over the station without drilling or modification to the station. This plate is NOT USED with the R5KLS-4.




Dimensions: Cover Plate: 4.49"W x 4.70"H


Back Plate: 4.25"W x 4.25"H

Back Box: single gang masonry style (will not fit ganged cut-in box), 2 ½ deep


Ronco High Security Solution for Responder 5 Nurse Call System

 RONCO <small>184 BRADY AVENUE HAWTHORNE, NY 10882</small>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 10%;">REV</th> <th style="width: 60%;">DESCRIPTION</th> <th style="width: 30%;">DATE</th> </tr> <tr> <td> </td> <td> </td> <td> </td> </tr> </table>	REV	DESCRIPTION	DATE				<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 10%;">DATE</th> <th style="width: 80%;">DESCRIPTION</th> </tr> <tr> <td> </td> <td> </td> </tr> </table>	DATE	DESCRIPTION			<p>RESPONDER 5 HIGH SECURITY SOLUTION HIGH SECURITY DRAWINGS</p> <p style="font-size: small;">/HAWTHORNE, NY</p>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;"> PART SYMBOL QTY UOM PART NAME </td> <td style="width: 80%;"> KKH NONE RS HIGH SECURITY 20.11.18.DWG </td> <td style="width: 10%;"> DATE 11/18/2020 INSTALLATION </td> </tr> </table>	PART SYMBOL QTY UOM PART NAME	KKH NONE RS HIGH SECURITY 20.11.18.DWG	DATE 11/18/2020 INSTALLATION
REV	DESCRIPTION	DATE															
DATE	DESCRIPTION																
PART SYMBOL QTY UOM PART NAME	KKH NONE RS HIGH SECURITY 20.11.18.DWG	DATE 11/18/2020 INSTALLATION															


STEP 1, SNAP STATION INTO BACK MOUNTING PLATE.




STEP 2, PLUG IN STATION, INSERT SCREWS THROUGH STATION AND BACK MOUNTING PLATE. SECURE INTO BACK BOX.




STEP 3, INSTALL PUSH BUTTON FACE PLATE.




STEP 4, ALIGN HIGH SECURITY COVER PLATE. ALIGN HIGH SECURITY COVER PLATE. SECURE TO STAND OFFS ON BACK PLATE.




STEP 1, SNAP STAFF STATION FACE PLATE ONTO STAFF STATION.



STEP 2, FROM BACK OF STAFF STATION, DRILL THROUGH FACE PLATE TO OPEN THE FOUR MOUNTING HOLES. (BOTTOM HOLES DO NOT NEED TO BE SLOTTED)



STEP 3, ALIGN HIGH SECURITY COVER PLATE OVER STATION. INSERT THE FOUR MOUNTING SCREWS THROUGH COVER, FACE PLATE AND STATION. SECURE INTO BACK BOX.



R5WPPB Vandal Proof/Water Resistant Call Switch



Vandal Proof/water resistant call switch is use in shower areas and wet locations that require a vandal proof switch. The switch housing and actuator are made of a high-quality stainless steel and is particularly suitable for use in harsh environments. Switch meets IP protection class IEC/DIN/EN/60529 (IP67). The switch is mounted on a stainless steel face plate and requires a one gang (Gangable or Non-gangable), 2 1/2" Deep back box. When installing, locate the switch away from direct water spray.