



**sepmag®**

## Qualitance v3.2

Biomagnetic Separation Monitoring  
& Management Software

### Quick Guide

Rev 20200505

**sepmag®**


## Qualitance

- The “Qualitnce” software monitors every single batch, allowing batch-tracing and the early detection of potential quality issues.
- It should be used in combination with the hardware is integrated in all Sepmag systems and its inserts



## Qualitance software

- Allows definition of multiple step processes.
- Identifies each **Sepmag** individually.



## Introduction

- **Qualitance** measures the changes on transmitted light through the bottle.
- At the start of the process, when the suspension is homogenous, the opacity is maximal.
- When separation is complete, the remaining suspension is clear, the opacity is minimal.



## Processes while the vessel is inside a Sepmag


During the time the vessel is inside the Sepmag:

- Biomagnetic separation takes place.
- The supernatant is extracted.
- Clean buffer is added to the vessel.

Every step should be performed when the previous one is completed.

**Qualitance records all the processes, but only displays Biomagnetic Separation.**

The complete record is accessible when the data is 'exported'.




## Start/stop



- **Qualitance** measures and records the changes from the introduction of the bottle inside the **Sepmag** until the extraction.
- The **Qualitance** is triggered automatically when the bottle enters the **Sepmag** and it does automatically stop when the vessel is removed.
- The 'End' or 'Start' can also be 'forced' by pressing the screen keys.

## Biomagnetic Separation Process

- The Biomagnetic Separation process occurs between the introduction of the bottle and ' $t_f$ '. The software ONLY displays this process (but records data until the bottle extraction)
- **Qualitance** software fits the experimental values to a sigmoidal curve between  $t=0$  and  $t=t_f$ , determining  $t_{50}$  and the exponent  $p$ .
- ' $t_f$ ' is determined by the user during the first measurement and can be edited through the Administrator menu.




May 5th, 2020

Qualitance v3,2

7

## $t_{50}$ and $p$ data management

- The software fits the experimental results and determines the Biomagnetic Separation parameters ( $t_{50}$  and  $p$ ).
- It also manages the historic data to generate reference curves for each Process




May 5th, 2020

Qualitance v3,2

8

## Predictive curve

- During measurement, once  $t_{50th}$  is reached, the software shows a predictive curve using the reference values ( $t_{50th}$ ,  $p_{th}$ )



## Content:


- Installation
- Quick Measurement
- Printing a Production Report
- Comparing different steps.
- Generating a new Process Macro
- First measurement of a new defined Process
- Selecting/Deselecting curves for generate the standard curve

## Installation:

- Download the installer
- Run the installer (\*.exe)
- Click on the 'Sepmag' button that would appear on the Desktop
- The default **username** is admin
- The default **password** is admin


## Quick Measurement (1/4)

- Make sure the **Sepmag** is connected to an USB port of the computer running the software.
- Log into the program.
- Click on Measures/'Process Monitoring'





## Quick Measurement (2/4)



- If the **Sepmag** is connected, **Ready** will be displayed on the screen.
- Select a process (a 'Standard Process' is included by default).
- Write the name of the Lot


May 5th, 2020

Qualitance v3,2

13



## Quick Measurement (3/4)




- Select the Step and the associated **Sepmag**.
- **Waiting**, will flash on the screen.
- The system is ready for inserting the bottle.

May 5th, 2020

Qualitance v3,2

14

## Quick Measurement (4/4)



- Put the bottle inside the **Sepmag** and the measurement will start automatically.
- When the measurement reaches the defined 'final time' a beep sound will play.


May 5th, 2020

Qualiance v3,2

15

## Printing a Report

- Go to Measure/Retrieve Data
- Select a Lot, then 'View Lot'
- The graph will show the experimental data and the reference curve. This last, can be hidden



May 5th, 2020

Qualiance v3,2

16

## Comparing different steps (1/2)



- When developing new process, it can be interesting compare the current step with the previous one, or with other similar process.
- Multi selection curve option, allows to compare several 'steps'.
- Go to 'Retrieve Data' screen.


May 5th , 2020

Qualityance v3,2

17

## Comparing different steps (2/2)

- Select the Product, Lot and Step, and click 'Add'.
- Click 'View' to 'show' the graphs (can be 'selected'/'unselected' at the table)
- Steps can be removed, selecting the row and clicking 'Remove'
- Print would generate a 'printer friendly report'



May 5th , 2020

Qualityance v3,2

18

## Generating a new Process Macro (1/2)

The screenshot shows the 'Qualityance' application window. The 'Administration' menu is open, highlighting 'Add/Modify Processes'. The 'Modify Processes' dialog box is displayed, showing a table of existing processes:

Product Name	Steps	Last User Modification	Date Last Modification
LAB	1	admin	2017/05/16 12:25
NEW PROCESS	2	admin	2017/05/16 13:22
Standard Process	1	admin	2017/01/02 12:39

Below the table, there is a form for adding a new process. The 'Product Name' field is set to 'CUSTOM PROCES' and the 'Steps' field is set to '2'. There are buttons for 'Update Product', 'New Product', 'New Step(s) for Selected Step', and 'Select Curve'. A 'Close' button is at the bottom right.

- Go to Administration/ Add/Modify Processes
- Introduce the name of the Product and the number of Steps

May 5th , 2020

Qualityance v3,2

19

## Generating a new Process Macro (2/2)

The screenshot shows the 'Add/Modify Processes' dialog box. The 'Step' table is visible, with columns for Step, Step Name, t50, P, if, QCR, Sampling period, and Trigger (s.u.). The 'Product Name' field is set to 'CUSTOM PROCES' and the 'Steps' field is set to '2'. There are buttons for 'Update Product', 'New Product', 'New Step(s) for Selected Step', and 'Select Curve'. A 'Close' button is at the bottom right.

Step	Step Name	t50	P	if	QCR	Sampling period	Trigger (s.u.)
1		0	0	0		100 ms	150
2		0	0	0		100 ms	150

- In the table, assign a name to the steps, and associate a SEPMAG device to each.
- Click on 'New Product'

May 5th , 2020


Qualityance v3,2

20



# First measurement of a new defined Process (1/3)

- When a new process is defined, the first measurement doesn't have the parameters defined.




May 5th , 2020

Qualitance v3,2

21



# First measurement of a new defined Process (2/3)




- It is important to click 'Force End' when the separation is complete. This will define the final time  $t_f$  to calculate the process parameters.

May 5th , 2020

Qualitance v3,2

22

## First measurement of a new defined Process (3/3)



- For the second measurement and the following,  $t_{50}$ ,  $p$  and  $t_f$  are already defined and a predictive curve will be drawn during the process (once  $t=t_{50}$ ).
- The system will be measuring until the vessel is removed from the Sepmag.


May 5th , 2020

Qualitance v3,2

23

## Selecting/Deselecting curves for generating the reference curve (1/2)

- By default, the software will calculate the 'theoretical'  $t_{50}$  and  $p$  for each process, taking the average of the values of all the experiments.
- The Administrator can remove curves (if some of the measurements were wrong) by clicking 'Select curves' at the Processes Macro screen.



May 5th , 2020


Qualitance v3,2

24

## Selecting/Deselecting curves for generating the reference curve (2/2)

- Curves can be add/removed just by clicking.
- Any change will be recorded in the Log File

- Once the statistics are good enough, the results can be 'frost'. Then no data from new curves will be added.



- For noisy data (low concentration suspensions), a running average filter may be activated.
- The number of initial/final excluded points can be also modified to improve the fitting quality

May 5th , 2020

Qualitance v3,2

25

## More information...

For additional information, questions and technical support, please contact your local sepmag® representative or email to [contact@sepmag.eu](mailto:contact@sepmag.eu)

For additional resources on biomagnetic separation (free eBooks, technical Posts, etc.), visit our website [www.sepmag.eu](http://www.sepmag.eu)

May 5th , 2020

Qualitance v3,2

26