

SARS-CoV-2 (N gene detection)
Test Information for Laboratory Manuals (United States)

Test name:	SARS-CoV-2 (N gene detection)
Test Description:	Qualitative real-time reverse transcription polymerase chain reaction
Specimen Collection / Specimen Type:	<p>Nasopharyngeal (NP) specimen is preferred. When collection of a nasopharyngeal swab is not possible, the following are acceptable alternatives:</p> <ul style="list-style-type: none"> • An oropharyngeal (OP) specimen collected by a healthcare professional, or • A nasal mid-turbinate (NMT) swab collected by a healthcare professional or by onsite self-collection (using a flocked tapered swab), or • An anterior nares specimen (NS) collected by a healthcare professional or by onsite self-collection (using a round foam swab). • For NS, a single polyester swab with a plastic shaft should be used to sample both nares. <p>Specimens should be collected into viral or universal transport media, Amies, or RNase free Normal saline.</p>
Acceptable Volume:	1mL (0.8mL minimum)
Days/Times of Analysis:	Sunday – Saturday, variable times
Specimen TAT:	<p>72 hours from receipt of specimen at Exact Sciences Laboratories (ES Labs)</p> <ul style="list-style-type: none"> • For positive results, ES Labs will place a call to the ordering provider at the phone number indicated on the test requisition form (as applicable). • ES Labs will report results to state health departments in accordance with communicable disease reporting requirements.
Transport:	Refrigerate specimens after collection and during transport to ES Labs
Stability:	<p>Specimens are stable at refrigerated temperatures for up to 72 hours after collection.</p> <p>Specimens can be stabilized at -70C for up to one month for clinical testing.</p> <p>Specimens should be shipped frozen on dry ice if they are anticipated to arrive at the lab beyond 72 hours from collection.</p>
CPT Code:	U0003
LOINC Code:	94533-7 SARS coronavirus 2 N gene [Presence] in Respiratory specimen by NAA with probe detection
McKesson Z-Code Identifier	None
Reference Value:	The normal value for this test is “Negative”

<p>Test Results:</p>	<p>Positive – SARS-CoV-2 (COVID-19) RNA detected.</p> <p>Negative –SARS-CoV-2 (COVID-19) RNA NOT detected. Negative results do not preclude SARS-CoV-2 (COVID-19) infection and should not be used as the sole basis for treatment or other patient management decisions</p> <p>Invalid – This specimen exhibited inhibition in the PCR assay or the specimen contained an inadequate amount of clinical material. Repeat testing is suggested if clinically indicated.</p> <p>Inconclusive- This specimen did not meet the full criteria established for the detection of SARS-CoV-2 (COVID-19) RNA. The sample will be evaluated using an alternative test method</p>
<p>Methodology:</p>	<p>Real-time RT-PCR test intended for the qualitative detection of nucleic acid from SARS-CoV-2 in respiratory specimens. Two primary steps:</p> <ul style="list-style-type: none"> • Extraction of viral RNA from patient specimens • One-step reverse transcription and PCR amplification with SARS-CoV-2 specific primers & real-time detection SARS-CoV-2 specific probes for the N gene of the SARS-CoV-2 viral RNA and the human RNase P (RP) gene.
<p>Analytical Sensitivity:</p>	<p>2.6 copies per µl specimen</p>
<p>Disclaimer Statement:</p>	<p>The SARS-CoV-2 (N gene detection) Test is a real-time RT-PCR assay intended for the qualitative detection of nucleic acids from SARS-CoV-2 in nasal, nasopharyngeal, and oropharyngeal swab specimens, from individuals suspected of COVID-19 by a healthcare professional. Testing is limited to Exact Sciences Laboratories, certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. § 263a, to perform high complexity tests.</p> <p>Results are for the detection of SARS-CoV-2 RNA. The SARS-CoV-2 RNA is generally detectable in respiratory specimens during the acute phase of infection. Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease</p> <p>Laboratories within the United States and its territories are required to report all positive results to the appropriate public health authorities.</p> <p>Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.</p> <p>This assay is intended for use under the Food and Drug Administration’s Emergency Use Authorization.</p>
<p>Aliases</p>	<p>COVID-19, nCoV-2019</p>