## Dräger Babyleo® TN500 IncuWarmer

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## Ventajas

---

### Redefining Thermoregulation – Keeping the temperature stable

Babyleo is Dräger's first IncuWarmer that provides optimal warming as an incubator, as a radiant warmer and during the transition between closed and open care. The convective heater, dual radiant warmer and heated mattress are synchronized to keep the temperature stable so heat loss and overheating are avoided. The convective heater is activated on a low level while the radiant heater is running so when you decide to switch modes, the transition can happen immediately giving you quick access without having to cope with heat loss. Advanced therapy warm-up, weaning and cooling modes allow you to focus on the patient. The closed humidification system allows you to set the humidity level up to 99% RH in manual mode or gives you the option to use Auto mode.

---

### Making your workflow easier and more efficient

The hood design allows you to open it with only one finger giving you exceptionally fast and easy access. The jumbo handports enable you to reach the entire bed area. The variable height adjustment ensures a comfortable workplace and allows mothers to be closer to their babies, even when sitting in a wheelchair. You can intuitively navigate through the display with the easy-to-use Dräger user interface and Touch-Turn-Confirm operating concept. Automatic modes are included so you have more time to focus on the baby instead of the device.

---

### The next best thing to the womb

Designed with the baby's development in mind, Babyleo has an extremely low operating noise level as well as light and noise monitoring inside the patient compartment. This way you are aware of the baby's exposure to potentially harmful stimuli and you can react to reduce these stressors. A safely integrated audio stimulation function gives you the opportunity to play mom's heartbeat or voice to soothe and stabilize the baby.

---

### Welcoming family centered care

The friendly device design, large variable height adjustment range plus other elements allow parents to comfortably get closer to the baby and make families feel welcome when visiting. The dedicated and personalized Family View clearly displays baby's therapy status and progress so it is easy to understand at a glance – all in an effort to help integrate parents into the care process.

---

### Designed with infection prevention in mind

The device's smooth surfaces and limited number of parts for disassembly make cleaning fast and easy. The completely closed humidification system, plus the fact that humidification is provided by steam from boiled water, eliminates the risk of germs entering the system. The automatic humidifier cleaning mode heats up the remaining water in the humidification reservoir until the reservoir is clean and empty when the device is not in use.

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Ambient Conditions During Operation	Value
Temperature	20 to 35 °C (68 to 95 °F)
Ambient pressure	620 to 1,100 hPa (9.0 to 16.0 psi)
Relative humidity	20 to 95 %, non-condensing
Power Supply	
Mains power connection	100 V to 240 V, 50/60 Hz
Maximum power consumption	1000 W
Dimensions	
Device (width x depth)	1,154 mm x 690 mm (45.43 in x 27.17 in)
Height (with hood closed)	1,850 to 2,250 mm (72.83 to 88.58 in)
Height adjustment range	400 mm (15.75 in) Electrical height adjustment for optimum adaptation to clinical users
Weight	<140 kg (308.6 lbs) (without options and accessories)
Mattress area	450 mm x 690 mm (17.72 in x 27.17 in)
Adjustment range mattress	90 mm (3.54 in)
Inclination of the mattress tray	Continuously variable electrical mattress tray tilt mechanism (up to 13°) to both sides and bed height adjustment
Side wall height	205 mm (8.07 in) with mattress tray in bottom position
X-ray tray	464 mm x 379 mm (18.27 in x 14.92 in) For larger film formats remove the x-ray tray and insert detector. Tray is accessible from outside for taking x-rays without disturbing the patient.
Display	264 mm (10.4 in) Color-Touchscreen
Operating Modes	Measured Value Display
Air temperature regulation	Setting range 20 to 39 °C (68 to 102.2 °F) <28 °C (82.4 °F), after confirmation >37 °C (98.6 °F), after confirmation Measurement range 13 to 45 °C (55.4 to 113 °F)
Radiant warmer	with 2 heating elements for an optimal and even heat distribution on the mattress area, automatic adaption of radiant warmer when mattress is tilted Setting range Off, 10 % to 100 % Pre-warm mode: 100 % for 3 minutes 60 % for 11.5 minutes 30 % until setting is changed
Skin temperature regulation	Setting range 34 to 38 °C (93.2 to 100.4 °F) >37 °C (98.6 °F), after confirmation 13 to 43 °C (55.4 to 109.4 °F)
Skin temperature monitoring	2 temperature sensors for central and peripheral skin temperature, graphical display of skin temperature
Operating Mode	Measured Value Display
Heated Mattress (Optional)	X-ray translucent, x-ray tray can be used for acquisition of x-ray images Setting range Off, 35 to 39 °C (95 to 102 °F), AUTO
Humidifier regulation	Closed servo-humidification system with external water reservoir, boiling and evaporation of Aquadest.

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Oxygen regulation (Optional)	Water bag holder on rear of main column. Normal measurement range: 30 to 99 % RH Setting range: Off, 30 to 99 % RH in steps of 1 %, AUTO Servo-controlled Oxygen enrichment in patient compartment with connection to central O <sub>2</sub> supply and/ or to O <sub>2</sub> -cylinder mounted on device
Flow velocity over mattress surface	Measurement range Normal range: 18 to 65 Vol% Extended range: 66 to 99 Vol% Measurement inaccuracy in normal range: ± (2,5 Vol% + 2,5 % from reading) O <sub>2</sub> rise time from 21 to 65 Vol% <10 minutes
Operating volume in patient compartment	10 cm/s (3.94 in/s) (with a horizontal mattress tray and an air temperature of 36 °C (96.8 °F))
Noise measurement in patient compartment (Optional)	Typically 40 dB(A), measured without oxygen application
Light measurement in patient compartment (Optional)	Measurement range 40 to 99 dB(A)
Scale (Optional)	Measurement range 3 to 999 Lux (1 to 93 fc) Measuring range (incl. tare): 200 to 10,000 g (0.44 lb to 22 lb) Tare weight: <= 10kg (22 lb) Display range: 0 to 10,090 g (0 to 22 lb 2 oz) Resolution (NAWI version): 10 g (0.022 lb) Resolution (standard version): 1 g (1 oz)
<b>Ergonomic Access to Patient</b>	
Hood	Detail Quick and simple mechanical opening and closing of the hood, power supply independant, also possible when device is switched off, softened closing, opening time typically < 2 sec
5 Hand ports	2 Jumbo-Handports per side for access to entire mattress area plus additional hand port at the front
2 side access panels and front access panel	Hinge down all access panels for patient access for medical interventions and primary care. Pulling out the bed support for medical procedures and for parents enabling close contact to their baby.
10 Grommets	Double at each corner and single grommets at front and back of the patient compartment
<b>Further Functionality</b>	
Views	Detail Main page, Thermomonitoring, Developmental Care, Family View
Graphical and tabular display of trends	2 hours to 7 days for temperatures, heating, noise level/ light level measurement, weight, events
Lock	Screen, height adjustment, bed tilt against accidentally adjusting the device settings
Timer	Stopwatch and APGAR timer with alarms after 1 min, 5 min, 10 min.
ClearView	Condensation management, the radiant warmer runs on a low level to reduce condensation on the inside of the hood.
Transition	Standby operation between incubator and radiant warmer to ensure a rapid change of operating mode without loss of heat. In incubator operation, the radiant warmer is running on a low level, in radiant warmer operation the convection heater is pre-heated.

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Alarm management	Ramping alarm volume, suppression and delaying of alarms automatically and proactively, Audio pause, AutoSet for alarm limits, Alarm volume setting range 50 – 70 dB(A)
Playing Audio files (Optional)	Connect external MP3 Player or cellular phone for playing music or parents' voices via the integrated loudspeakers, volume is limited to a maximum of 55 dB(A) regardless of the volume set.
Light	Working light (3 intensity options), night light, moodlight (7 color choices available).
Storage (Optional)	Lockable drawer with insert
Castors	Arms for monitors, IV pumps, examination lights and shelves can be mounted on the duct rails on right and left side of the main column. Standard rails on front and rear for mounting accessories and bed coupling.
Reprocessing	4 double castors with locking brakes, one can be locked as directional castor for straight movement
Data Exchange and Interfaces	Simple device reprocessing, only a few number of parts have to be disassembled, plain surfaces, surface disinfection of all parts Automatic humidifier cleaning Serial port COM (Medibus.X Protocol), USB port, Nurse call, Service port (RJ45)
<b>Materials</b>	<b>Detail</b>
Bisphenol A (BPA)	The patient compartment of the device (without accessories) is BPA-free
Latex	The device (without accessories) contains no latex-containing components
Gas cylinder holder (Optional)	Quick release holder for attaching gas cylinders (1 oxygen and 1 air) to the trolley <ul style="list-style-type: none"> <li>– max. total height 33 – 76 cm (12.99 to 29.92 in)</li> <li>– max. diameter 10 – 14 cm (3.91 to 5.51 in)</li> </ul>

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