



UrgoKTWO

LEG ULCERS



UrgoKTWO how to apply

For venous and mixed aetiology leg ulcers

Below the knee application using 10cm width bandage kits

UrgoKTWO KITS CONSISTS OF: LAYER 1: wadded short stretch bandage **LAYER 2:** Long stretch cohesive bandage (latex free available)

Before bandaging, measure the ankle and choose the correct **UrgoKTWO** kit (between 18-25 cm or 25-32cm). **UrgoKTWO** (~40mmHg*) for venous ulcers and **UrgoKTWO Reduced** (~20mmHg*) for mixed aetiology.

A. BANDAGE APPLICATION

- Position foot at 90° angle (**toes to nose**)
- Ensure the wadding side of layer 1 is next to the skin. Anchor in place with 1-2 turns at the base of the toes
- Secure the heel with a figure of 8. Ensure foot and heel are fully covered. Continue bandaging in a spiral up the leg. **Stretch the pressure indicator to form a circle**
- Overlap the bandage so the pressure indicator is just covered at their base (see section B)
- Finish bandaging 2cm below the Popliteal space
- Cut off excess bandage and secure with tape
- Apply layer 2 using the same technique
- For patient comfort, allow a small border of layer 1 at the toes and knee. Cut off any excess and press down gently to ensure cohesion

B. POSITION OF PRESSURE INDICATORS

 Size: 18-25cm Position of pressure indicator: Top edge 50% overlap	 Size: 25-32cm Position of pressure indicator: Middle of bandage 2/3 overlap
--	---



UrgoKTWO Reduced / UrgoKTWO Reduced Latex Free (~20mmHg*) has a R in the pressure indicator to differentiate from UrgoKTWO (~40mmHg*)

*average donated pressure

Carefully read the leaflet before use.



Validates the TLC Healing Matrix

191150 August 2018