

# Intrafix® Primeline Classic

The standard IV set



Meets all your routine fluid administration needs with state-of-the-art technology and great value

**B | BRAUN**  
SHARING EXPERTISE

# Meet today's challenges ...

## Intrafix® Primeline Classic – ergonomic in design, enduring in tradition

The universal practicality of this classic in the Intrafix® line has been proven billions of times worldwide. User-friendly core elements like the drip chamber and roller clamp enable secure gripping and easy handling.

B. Braun's high-end production technology and unrivaled degree of automation guarantee consistently dependable high quality.

Intrafix® Primeline Classic offers everything a standard IV administration set needs – perfect to the last detail.



Clear, strong ISO spike

Highly transparent dome design

Unique center ring

Elastic pump chamber

Large filter area

See-through tubing  
150 cm in length

Screw-tight  
lock fitting

Handy roller clamp

# ... the well-established way

## Fluid Administration

### User benefits

- Secure gripping
- Effortless and safe spiking
- Swift setting
- Easy monitoring
- One-handed adjustment

### One spike fits all

The universal, vented spike fits tightly on all container types. Removal is just as smooth as insertion. Fluid inlet and air outlet are separated at different levels to prevent inflowing air from inadvertently entering the fluid channel.



### Instant drop speed adjustment

The large-dimensioned roller clamp enables secure gripping. The spike holder is an added feature that protectively covers the tip after use. The tubing can be hooked into the side notch after priming for stand-by procedures.



### Universal and correct fit

Insert the connector tip, screw in tightly with a quick twist. Prevents inadvertent disconnection.

### Solve the "spaghetti syndrome"

No tangling, thanks to the tubing's perfect length. Minimizes residual infusion fluid loss.

### Ergonomically designed

Larger and longer than conventional, this drip chamber is ideally shaped to fit firmly in your hand. Its center ring ensures secure gripping. The soft, flexible lower pump chamber enables one-click fluid level setting.

The fully transparent, dome-shaped upper part allows a quick drop rate check even in poor lighting.

The large filter membrane retains  $>15\ \mu\text{m}$  particles.

Pump drop detectors are easily attachable.



# Choose your best Intrafix® option

## Fluid Administration

Intrafix®	Specification	Length of tubing	Code No. (REF)
SafeSet	Basic	180	4063000
	B.C.V. (with Back Check Valve)	180	4063001
	N.T.P. (with Neutrapur, PVC-free)	180	4063002
	Basic	230	4063003
	Y-N.F. (with Y-site Needle-Free Valve)	210	4063004
	I.S. (with Injection Site, with Grip Plate)	180	4063005
	3-W.S.C. (with 3-Way Stopcock)	220	4063006
Primeline Comfort	Basic	180	4062981L
	PrimeStop	180	4062982L
	B.C.V. (with Back Check Valve)	180	4063287
	Basic	230	4060407
	N.T.P. (with Neutrapur®, PVC-free)	180	4062191
	I.S. (with Injection Site, without Grip Plate)	180	4062181
	I.S. (with Injection Site, with Grip Plate)	180	4062182
	Y-N.F. (with Y-site Needle-Free Valve)	180	4062158
	Air-matic®	180	4063988
	Air-matic® B.C.V. (with Back Check Valve)	180	4063783
Primeline Classic	Basic	150	4062957

Sales unit / Pcs.: 100



**B | BRAUN**  
SHARING EXPERTISE

**B. Braun Melsungen AG**

34209 Melsungen  
Germany  
Tel +49(0) 56 61 71-0

[www.bbraun.com](http://www.bbraun.com)