

3point products® 3pp® Design Line™ Thumb Splint



The **3PP DESIGN LINE THUMB SPLINT** combines comfort and fashion in a whole new way. Designed to support the CMC joint for relief of the grinding pain of CMC joint arthritis, it's easy-to-apply and fully adjustable for a customized fit. Available in 4 fashion patterns, this moderate support splint provides **COMFORT, FIT** and **FASHION**.

Features

- Contoured strap supports and carefully abducts the CMC joint to reduce the pain from CMC joint arthritis or ligament injury.
- Soft foam-lined stretch material contours to the hand with natural compression so there is no need to overtighten the splint.
- Straps can be secured anywhere on the hook-engageable material to achieve a custom fit.
- Three sizes and three patterns offer an ideal range of fit and fashion and are the perfect alternative to bulky, black-on-black splints.
- 100% breathable, 100% latex-free; hand or machine washable.

Order Information

Measure across the palm at the level of the big knuckles (MP joints).



Size
Measure across the palm at the level of the big knuckles (MP joints)

Small 2 3/8" to 3" (6cm to 7.6cm)
Medium 3" to 3 5/8" (7.6cm to 9.2cm)
Large 3 5/8" to 4 1/2" (9.2cm to 11.4cm)

Hand & Size	Order #	Style Choice
Left Small	P2008-L2	Black Camo add BBC
Left Medium	P2008-L3	Red Trim add RB
Left Large	P2008-L4	
Right Small	P2008-R2	
Right Medium	P2008-R3	
Right Large	P2008-R4	

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com