




## Helping Hands Heal for Over 25 Years





## Splints Designed for Comfort, Fit and Function



Original Oval-8® Finger Splints help manage Mallet Fingers, Swan Neck, Trigger Finger and Thumb, Boutonniere Deformity, Hypermobility and more.

Designed with 3-points of control for optimal joint support. Oval-8s are effective, comfortable, and easy to wear day or night.

The patented band allows each of the 14 sizes to fit tighter or looser depending on which end is put on first.

Available in CLASSIC BEIGE and CLEAR



Boutonniere



Swan Neck



Trigger Thumb



Mallet



### 3pp® Buddy Loops®

Fast, Comfortable Control



Simple loop and wrap design and non-slip foam lining make 3pp Buddy Loops the best option to protect and restore motion.



Oval-8® Finger Splints, 3pp® Buddy Loops® and 3pp Toe Loops® are registered in Eudamed and approved for EU distribution.

No natural rubber latex is used in any 3-Point product or packaging.

3 widths:  
1/2" (1.27 cm)  
3/4" (1.9 cm)  
1" (2.5 cm)



## FINGER SPLINTS



### 3pp® Finger Trapper



Provides comfortable, non-slip traction for splinting or post-op positioning. Conforms over the nailbed to eliminate pressure. **One size fits most.**

### 3pp® Final Flexion Wrap



Restores end range PIP & DIP flexion. Easy to apply static progressive stretch. **One size fits most.**

### 3pp® Side Step Splint



Adjustable stretch corrects IP joint deviation caused by arthritis or injury. Adjust tension to tolerance to improve alignment. **3 sizes. Not for use on thumb.**

### 3pp® Step Down Splint



Restore composite DIP and PIP flexion following injury or surgery. Easy to apply tension to tolerance for optimal results. **3 sizes. Not for use on thumb.**

## ULNAR DEVIATION SPLINTS

### Polycentric Hinged Ulnar Deviation Splint



Holds MP joints in neutral and allows full flexion and extension. Recommended for deviation that's correctable with mild to moderate force. **4 sizes in right and left models.**

### Radial Hinged Ulnar Deviation Splint



Each finger is adjusted individually to improve alignment of strongly deviated fingers. Radial hinge allows free flexion and extension. **3 sizes in right and left models.**



### Comforter Splint



Softly padded splint holds the wrist and MP joints in neutral to relieve stress. Reduces morning stiffness for better motion and function. **3 sizes in right & left models.**



## THUMB SPLINTS

### 3pp® Prima Thumb Brace



Optimal support to treat CMC and MCP Arthritis. Features fully adjustable strapping, memory stretch material and removable stay. [3 sizes, right and left models.](#)

### CMCcare Thumb Brace



Long lasting, supportive thumb brace stabilizes the CMC joint for pain-free motion. Includes a cushioned liner, removable and replaceable straps. [3 sizes, right and left models.](#)



### Fix Comfort Thumb Brace



Lightweight, adjustable brace provides firm support for Ligament Injury or Arthritis. Fully cushioned in tri-layer foam, easy to wear all day and/or night. [2 sizes fit both right and left.](#)

### 3pp® Ez FIT Thumb Spica Splint



Fully cushioned in memory foam to protect de Quervain's and Gamekeeper's Thumb. Full length malleable aluminum stay provides firm support. [One size fits right or left hand.](#)

### 3pp® ThumSaver MP



Protects the thumb MP and stabilizes CMC motion without restricting wrist motion. [4 sizes, right and left models.](#)

### 3pp® ThumSlings



3pp ThumSling provides light support to control CMC pain. Choose the 3pp ThumSling NP for a greater degree of support. [4 sizes, right and left models.](#)

## WRIST & ELBOW SPLINTS

### 3pp<sup>®</sup> Carpal Lift



Helps manage Ulnar Carpal instability with a unique strap and pad system. Relieves pressure on Guyon's Canal without limiting wrist motion. [2 sizes, right and left models.](#)

### Fix Comfort Wrist Brace



One-handed strap system creates custom fit for wrist immobilization to treat Carpal Tunnel, wrist instability, ulnar impact, post-op, or cast removal. [Two sizes fit both right and left hand.](#)



### 3pp<sup>®</sup> Choice Wrist Brace



One-size fits most brace for Carpal Tunnel, impact injuries, and post-op or post-cast use. The dorsal stay can be adjusted to fit a wide range of sizes. [1 size, right and left models.](#)

### 3pp<sup>®</sup> Wrist P.O.P. Splint



Cushioned pads apply counterforce pressure on the distal radius and ulna to stabilize the DRUJ. Easy to adjust. [2 sizes, right and left.](#)

### 3pp<sup>®</sup> Elbow P.O.P. Splint



Adjustable pressure pad reduces pain from medial or lateral tendinosis. Non-slip foam lining holds the strap in place. [2 sizes, right or left.](#)

### Gel Mate Silicone Gel Sheeting



Medical grade silicone sheets are self-adherent, reusable and long lasting. Available in Smooth Beige and CLEAR. [4" x 6" sheets.](#)



### Gel Tubes



Protect and cushion finger lacerations, enlarged joints and more. [Use on fingers or toes.](#)



# Where to Purchase

Products are available through your contract suppliers as well as through major websites. Health Care Providers can buy direct from 3-Point Products at professional.

## Health Care Providers



888-378-7763  
Mon-Fri 8:30am -4:30pm EST  
service@3pointproducts.com



[www.3pointproducts.com](http://www.3pointproducts.com)

## U.S. Consumers and Patients

Our e-commerce store is an easy solution for your patients to purchase the products you recommend



[www.ShopOMA.com](http://www.ShopOMA.com)



## Stay informed. Stay updated. Stay connected.

3pp Clinic Connect offers health care providers live product demonstrations to learn about our products, our company, and how we can help your patients.

Sign up and attend our monthly live product demonstrations to learn about our products, ask questions, and receive FREE clinic samples.

Demos for Clinics & Group Practices  
Explore our products through customized demonstrations designed specifically for your practice.



Photo reproduced by kind permission of University of Utah Hand Therapy Clinic

