



3-Point Products Finger Orthoses

Product Name	Primary Indications	Splint Features	Availability
Oval-8 Finger Splint	Stabilize, immobilize or limit motion: <ul style="list-style-type: none"> • Mallet Finger • Boutonniere Deformity • Lateral Deviation • Fractures • Swan Neck Deformity • Hypermobility (Ehlers-Danlos) • Trigger Finger & Thumb 	<ul style="list-style-type: none"> • Long lasting molded polypropylene won't deform in normal heat • Splints fit tighter or looser depending on which end is slipped first • Use to immobilize, limit or restrict range motion for rest and healing • Heat adjustable for custom fitting; can be worn in water 	<p>Available in 14 sizes. Sold in single size and mixed size packages</p> <p>Oval-8 Kit, Sizing Set and Pediatric Kit are ideal for clinic sizing and dispensing</p>
3pp Side Step Splint	Reduce IP joint lateral deviation due to: <ul style="list-style-type: none"> • Osteo or Rheumatoid Arthritis • Collateral ligament injury • Heberden's or Bouchard's nodes 	<ul style="list-style-type: none"> • Static progressive splint lets wearer to adjust corrective force to tolerance • Adjustable, cushioned counterforce pad • Easy one-handed application • Steel wire frame will not deform with force 	<p>3 sizes - Small, Medium, Large</p> <p>Measure length and width and refer to sizing chart</p>
3pp Step Down Splint	Restore DIP & PIP joint flexion <ul style="list-style-type: none"> • Extensor tendon contracture • Capsular or lateral band adhesion • Interosseous muscle contracture 	<ul style="list-style-type: none"> • Static progressive splint lets wearer adjust corrective force to tolerance • Pads apply force to optimize PIP & DIP ROM • Easy one-handed application • Steel wire frame will not deform with force 	<p>3 sizes - Small, Medium, Large</p> <p>Measure length and width and refer to sizing chart</p>
3pp Buddy Loops	Protection and ROM assist for: <ul style="list-style-type: none"> • Jammed or sprained finger • Protect and assist range of motion • Post-op positioning • Relative motion strap 	<ul style="list-style-type: none"> • Loop and wrap-around design with hook closure • Easy one-handed application • Non-slip foam lining keeps loops in place • Three widths for a range of finger sizes & applications 	<p>Three widths - 1/2" (1.27 cm) 3/4" (1.9 cm) 1" (2.5 cm)</p> <p>All sizes are 5" (12.5cm) long</p> <p>Packages of 5, 25, 50 or 100</p>
3pp Final Flexion Wrap	Restore end range DIP & PIP flexion <ul style="list-style-type: none"> • Post extensor tendon repair • Capsular or lateral band adhesion • Interosseous muscle contracture 	<ul style="list-style-type: none"> • Easy to adjust static-progressive wrap-around splint • Non-slip foam lined fabric contours to the finger for comfortable pressure distribution • Ideal following use of Step Down Splint for final ROM 	<p>One size fits most</p> <p>Packages of 5 or 25</p>
3pp Finger Trapper	Application of finger traction for post-op dynamic splinting <ul style="list-style-type: none"> • Tendon repair splints • Dynamic or static progressive splinting 	<ul style="list-style-type: none"> • Non-slip foam lined straps fastens around the proximal • Eyelet for attachment of a static or dynamic traction line • Replaces fingernail sutures, hooks, or glue 	<p>One size (1/2" wide by 1 5/8" long) stretches to fit a wide range of finger sizes</p> <p>Package of 5</p>
Gel Tubes	Protect and cushion fingers affected by: <ul style="list-style-type: none"> • Finger lacerations or ulcerations • Cushion under finger splints • Helps soften and lessen scarring 	<ul style="list-style-type: none"> • Hypoallergenic mineral oil based polymer • Reusable, durable and washable • Latex-free 	<p>Available in single size and mixed size packages</p> <p>Oval-8 Kit, Sizing Set and Pediatric Kit are ideal for clinic sizing and dispensing</p>