



### 3-Point Products Thumb Orthoses Options and Features

Product Name	Indications	Splint Features	Support Level	Support Level Description
<b>3pp ThumSling</b>	CMC thumb arthritis Mild ligament sprains Overuse Injuries	<ul style="list-style-type: none"> <li>• Lightweight, slim-fit wrap around style</li> <li>• Soft non-slip foam lined material</li> <li>• Long contoured strap repositions and supports the CMC joint and wrist</li> </ul>	1	1= Light support uses compression & targeted support to control motion
<b>3pp Prima Thumb Brace</b>	CMC and MCP joint arthritis Sprains, overuse injuries Hypermobility (EDS)	<ul style="list-style-type: none"> <li>• Easy-on-Tab for one handed application</li> <li>• Memory-stretch fabric for better fit</li> <li>• Non-slip straps for custom sizing and support</li> <li>• Removable dorsal stay</li> </ul>	1	2= Moderate support limits motion to a safe range and supports joints to reduce stress
<b>3pp Design Line Thumb Splint</b>	CMC and MCP joint arthritis Sprains, overuse injuries Hypermobility (EDS)	<ul style="list-style-type: none"> <li>• Contours for a non-bulky support</li> <li>• Foam-lined fabric and strap target CMC support</li> <li>• Choice of colored binding or fabric design</li> </ul>	2	3 = Firm support protects targeted joint(s) but allows necessary hand function
<b>3pp ThumSling NP</b>	CMC and MCP thumb arthritis Mild sprains Overuses injuries Hypermobility thumbs(EDS)	<ul style="list-style-type: none"> <li>• Similar to ThumSling with more supportive materials</li> <li>• Long strap supports the CMC and wrist</li> <li>• Removable stay can be added as needed</li> </ul>	2	4 = Limits or immobilizes targeted joints for rest and protection to allow healing
<b>3pp ThumSling Long NP</b>	CMC Thumb Arthritis Ligament injuries Mild wrist sprains or instability	<ul style="list-style-type: none"> <li>• Combines adjustable wrist and thumb support</li> <li>• Removable stay can be positioned as needed</li> <li>• Stretch-Memory material for a secure fit</li> </ul>	2	
<b>Fix Comfort Thumb Brace</b>	Thumb Arthritis Gamekeepers or Skiers Thumb Overuse Injuries	<ul style="list-style-type: none"> <li>• Ultra-light thumb and wrist support</li> <li>• Adjustable frame applies firm CMC and wrist support</li> <li>• Straps attach anywhere for a custom fit</li> </ul>	3	
<b>CMCcare Thumb Brace</b>	CMC thumb arthritis Instability due to injury or reduced strength	<ul style="list-style-type: none"> <li>• Enhances CMC stability with a conformable orthosis &amp; malleable stay</li> <li>• Open palm design for full finger and thumb motion</li> <li>• Straps are adjustable and replaceable</li> <li>• Includes a cushioned 3pp Cozy Liner &amp; an extra strap and closure</li> </ul>	3	
<b>3pp ThumSaver MP</b>	Thumb arthritis Gamekeepers or skiers thumb Ligament injuries	<ul style="list-style-type: none"> <li>• Molded brace supports the CMC and MCP thumb joints</li> <li>• Felt pads included to cushion the plastic shell for comfort</li> <li>• Palm is free for hand function</li> </ul>	4	
<b>3pp Ez FIT Thumb Spica Splint</b>	De Quervains tenosynovitis Gamekeepers or skiers thumb Overuse injuries	<ul style="list-style-type: none"> <li>• Full length adjustable aluminum stay holds thumb securely</li> <li>• Angled straps lie flat on the arm for comfort and a secure fit</li> <li>• Wrist motion is limited but fingers are free for function</li> </ul>	4	