





www.3pp-eifu.com
English (EN)

How to Wear Oval-8® Finger Splints


1. Boutonniere Deformity



4. Mallet Finger




2. Crooked Finger




5. Swan Neck Deformity



3. Trigger Thumb





6. Trigger Finger



Two sizes in each splint

Oval-8 splints fit tighter or looser depending on which end is put on the finger first.



Adjusting an Oval-8

Oval-8 splints can be adjusted by a healthcare provider using a heating tool to change the angle, length and fit.



Distributors

UNITED KINGDOM

Performance Health UK
www.performancehealth.co.uk

Promedics
www.promedics.co.uk

AUSTRALIA

Performance Health
www.performhealth.com.au

OPC Health
www.opchealth.com.au


NEW ZEALAND

AT Therapy Limited
www.attherapy.co.nz



3-Point Products, Inc. U.S.A.
www.3pointproducts.com



3-Point Products, Inc.
118 Log Canoe Circle
Stevensville, MD 21666 USA
www.3pointproducts.com
+001.410.604.6393 




















Sallmann Medical GmbH
Alois-Stockinger-Str.16
5020 Salzburg Austria
service@sallmann-medical.eu



Swiss AR Services GmbH
Industriestrasse 47
CH-6300 Zug

3-Point Products' Symbols Legend EN

Symbols	Symbol Title	Symbols	Symbol Title
	Medical Device		Authorized European Representative
	CE Mark		Authorized Swiss Representative
	Consult instructions for use		Manufacturer
	Single patient multiple use		Date of Manufacture
	Recyclable Polypropylene Plastic		Country of Manufacture US
	Hand wash with water, diluted bleach or alcohol		Lot number
	Age Restriction. Not for children under 5. May cause choking if swallowed		Distributor
	Model number		Patient Information Website
	Latex free. Not made with natural rubber latex.		