



Oval-8® Finger Splint

Oval-8® Finger Splints stabilize and align the PIP and DIP Joints. This preferred finger orthosis treats mallet finger, swan neck and boutonniere deformities, trigger finger and trigger thumb, lateral deviation, fractures, arthritis, and hypermobility (EDS). Oval-8 Finger Splints are packaged for convenient use in any size clinic or treatment setting.



Oval-8 Kit



Oval-8 Sizing Set



Single & Mixed Sizes

Features

- Molded plastic splints have broad bands and rounded edges for comfort
- Angled band allows for tighter or looser fit to accommodate mild edema or changes in finger size due to swelling, temperature change or time of day
- Open-air design keeps skin clean and dry
- Can be worn in water when washing hands, showering or swimming
- Quickly and easily adjustable with a heat gun

Ordering Information

Packages of 1

End numbers correspond to splint sizes.

P1008-1-02	P1008-1-09
P1008-1-03	P1008-1-10
P1008-1-04	P1008-1-11
P1008-1-05	P1008-1-12
P1008-1-06	P1008-1-13
P1008-1-07	P1008-1-14
P1008-1-08	P1008-1-15

Packages of 5

P1008-5-02	P1008-5-09
P1008-5-03	P1008-5-10
P1008-5-04	P1008-5-11
P1008-5-05	P1008-5-12
P1008-5-06	P1008-5-13
P1008-5-07	P1008-5-14
P1008-5-08	P1008-5-15

Oval-8 Graduated Sets

RP1008-001	Sizes 2, 3, 4
RP1008-002	Sizes 4, 5, 6
RP1008-003	Sizes 6, 7, 8
RP1008-004	Sizes 8, 9, 10
RP1008-005	Sizes 10, 11, 12
RP1008-006	Sizes 13, 14, 15

Oval-8 Combo Packs

P1008-C1	Sizes 3, 4, 5, 6, 7
P1008-C2	Sizes 6, 7, 8, 9, 10
P1008-C3	Sizes 10, 11, 12, 13, 14

Oval-8 Kit	P1008-K
Oval-8 Pediatric Kit	P1008-PK
Oval-8 Sizing Set	P1008-S

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com