

## Sling safety checks to complete before every use

- ✓ Check for an in-date service tag
- ✓ Check the label is legible
- ✓ Check the stitching is intact
- ✓ Check for signs of wear and tear
- ✓ Check the straps and belts are working 100%
- ✓ Check for signs of excessive heat or stress

If you are unsure about the condition of your sling, you should take them out of service and contact your service and maintenance provider.

If you identify a fault with the sling or if the service is out of date, immediately report it to your service and maintenance provider and do not use the equipment.

This is not a substitute for manual handling training and is intended as a guide only