



If this is the goal

...this is the waterline treatment.

For your office, the patient comes first. That's why along with *simple and efficient*, you need an **effective** waterline treatment backed by science and proven in practice.

BluTube® is specifically engineered for greater effectiveness, which means improved water test results and more consistent CDC compliance.

THE
BRAND
NEW

BluTube®
Dental Unit Water Purification Cartridge



THEIR PROTECTION
YOUR ASSURANCE

Questions?

888.843.3343 | ProEdgeDental.com/BluTube



BluTube®

Dental Unit Water Purification Cartridge




Two is Better Than One

With over 30% of treated waterlines failing to meet the CDC standard for safe water (≤ 500 CFU/mL)[†], you need not just an efficient solution, but an effective one with a proven protocol.

BluTube® was designed with the latest research in treatment efficacy and specifically engineered for safer patients and better compliance testing results.

- › Two-cartridge system treats waterlines with 40% greater power.*
- › Safe for patient contact, non-toxic, non-allergenic, and non-corrosive.
- › Includes 6-month automatic replacement indicator for convenient reminder & "dummy straw" for easy shocking.
- › Installs in minutes & can be used with tap or distilled water.

Total Active Ingredient Release Over 1 Year



Item #	BluTube®	MSRP
BT180	Two 6-Month Cartridge System	\$279. ⁹⁵



THEIR PROTECTION
YOUR ASSURANCE

Order Today!

888.843.3343 | ProEdgeDental.com/BluTube

[†]Study cited in "Treating & Monitoring Dental Water", Compendium of Continuing Education in Dentistry, by Nancy Dewhirst, RDH, BS & John A. Molinari, PhD. Data collected from 22,196 consecutive waterline tests (R2A, heterotrophic plate count tests) conducted in 2017 by ProEdge Dental Water Labs. All product types provided by dental offices. Averages may not reflect your test results. All Pass/Fail measured against CDC's ≤ 500 CFU/ml standard. *Pass rate based on 40,947 consecutive waterline tests conducted in 2018 and 2019 by ProEdge Dental Water Labs. *Based on total active ingredient release over the course of 365 days or 240L following instructions for use (replacing BluTube cartridge at 6 months). The 2-Cartridge BluTube® system continuously treats for 365 days or 240L if water usage records are kept. Each individual cartridge should be used for 6 months or 120L if water usage records are kept.

©ProEdge Dental 7.2020_V3