


# Additional Vial Instructions

Additional vials may be added to any ProEdge Dental Water Labs testing kits.


## **1. Label Vials**

With a permanent marker, number each additional vial sequentially, starting with the next available number. For example, if you ordered a 6-Vial Test Kit, label your additional vial "7" and then "8" and so forth.



## **2. Collect Samples**

Fill sample vial to 90% capacity. Close lid securely.



**FLIP FOR STEP 3**

### 3. Complete Submission Form

Fill out the form **completely** to ensure accurate results. The vial number you put on the sample will be indicated on your results.



# Questions?

888.843.3343

ProEdgeDental.com

