Continuous Measure Results

- Click on the *reporting icon* in ENCOR Hospital Abstracted Measures
- From the drop-down menu Choose Data Quality
- Select Continuous Measure Results



- Set up Dataset Filters, enter the Date Range, choose selections & click Go.
- Select Hospitals for Multi-Hospital Facilities
- Select Individual Measure Set or Several Measure Sets.
- If you want to look at specific measures, you can search and select the measures by checking the box next to the Measures you want to include.
- Check the box to include *In-progress* cases.

Con	Continuous Measure Results									
۲	Dataset Filters									
F	rom Date	To Date		Hospitals	۲	Measure Sets	۲	Measures	•	Go
		h5 —		🗹 All		Select All		Select All		
				General 1		OP-SURGERY - OP-SUR	GERY - OP Surgery			
				Medisolv1		HRIDGLEUR,1 . HRIDGL	CUR.1 . HRIDG.CUR.1			
								Search		
1	Include In-progress Cases									

• By default, the report is grouped by Measure ID. The Mean and Median for each Measure is displayed.

Continuous Measure Results								
Dataset Filters								
Report Results								
Field List Expand/Collapse Excel PDF RTF Print Clear Filter								
L	Measure ID 🛆							
	Measure Set ID Measure	Measure Value	Patient Identifier	Abstract				
	V V	4	Ÿ	Ÿ				
15	Measure ID; ED-1a		Mean:370 Median:0					
E	i Measl्रुरोD: ED-1c		Mean:1,480 Median:1480					
E	Measure ID: ED-2a		Mean:353 Median:0					
E	Measure ID: ED-2c		Mean:1 411					
			Median:1411					
E	Measure ID: OP-18a		Mean:930 Median:930					
E	Measure ID: OP-18b							
			Mean:420 Median:420					
Ε	Measure ID: OP-18d		Mean:1,440 Median:1440					
1		Save View Layout Re	nove Saved View Layout					

• To view individual cases, click on the (+) symbol next to the Measure ID.

	Measure ID: OP-18a							
	OP-ED	Median Time from ED Arrival to ED Departure for Discharged ED Patients - Overall Rate	1537		2	4979		
	OP-ED	Median Time from ED Arrival to ED Departure for Discharged ED Patients - Overall Rate	2775		2	4980		
	OP-ED	Median Time from ED Arrival to ED Departure for Discharged ED Patients - Overall Rate	87	T.	2	4989		
	OP-ED	Median Time from ED Arrival to ED Departure for Discharged ED Patients - Overall Rate	420	L. L.	2	4991		

Identifying extreme values

- Identifying extreme values (those values that are larger than expected) can be accomplished by grouping by Measure Value and sorting Measure Value in descending order.
- Group by *Measure Value*, then sort Measure Value in ↓ Descending order by clicking the *Measure Value* label once, then a second time.
 - This will now list any potential extreme values at the top of the Measure Value column. An extreme value can be the result of an incorrect date or time entered.

Con	ontinuous Measure Results								
۰ ی	Dataset Filters								
۰	Report Results								
Field List Expand/Collapse Excel PDF RTF Print Clear Filter									
[Measure Value v Measure ID A								
	Measure Set ID Measure Patient Ident	tifier Abstract	Medisolv Identifier	Case Identifier					
	Ý	\$	· · · · · · · · · · · · · · · · · · ·	▼					
E	Measure Value: 11680								
Ē	Measure Value: 11491								
Ē	Measure Value: 1440								
Ē	Measure Value: 420								
1	Sav	e View Layout Remove Saved View	v Layout						

- To review the Extreme Values, drill-down on any Measure Value to view the Measure ID and Medisolv Identifier by clicking on the (+) symbol next to Measure Value.
- Click on the (+) symbol next to the Measure ID to expand your view.
- The abstract icon with standard color coding is displayed for each record and can be opened from this page.
 - Right click the abstraction icon to open in a new tab or new window or hold down the CTRL key and click the abstraction icon.

Measure Value ∇ Measure ID \triangle							
		Measure Set ID	Measure	Abstract	Medisolv Identifier		
		Ŷ	٩		v		
Measure Value: 41756							
		OP-CHEST PAIN	Median Time to ECG - Overall Rate		1719		

• This report can be set-up in the desired format using grouping and sorting. Once set-up, Save View Layout will make these selections the default layout the page opens in.



• The below pop-up appears after a View is saved.



• To remove a saved layout, click *Remove Saved Layout*. The below pop-up is displayed. Click *OK* to remove the saved view.

11/13/2020 - November	HRIPS-GLORAL	Swan
		×
This site says Are you sure you want to remove	e the saved view?	
ОК	Cancel	