



Qnostics

Respiratory Triplex Q Control

Respiratory Triplex Q Control

As a world leading manufacturer of third party controls for molecular infectious disease testing, the Qnostics range comprises hundreds of characterised viral, bacterial and fungal targets. Whole pathogen and manufactured independently, the Respiratory Triplex Q Control is designed to monitor the daily performance of molecular methods used in the detection of Influenza A, Influenza B and Respiratory Syncytial Virus A (RSV A).

Five positive and five negative samples are conveniently supplied in a single pack reducing the number of individual kits required. The target concentration of each of the positive samples is designed to be within the dynamic range of most molecular assays.

Key Benefits

- Whole pathogen control designed to monitor the entire testing process
- Liquid frozen for user convenience and ease of use
- Consolidated triplex control delivering significant cost and time savings
- Suitable for use with a wide range of commercial and in-house molecular methods
- Manufactured to ISO 13485 quality standards



Specifications

Target Pathogen – Influenza A (H1N1), Influenza B (Victoria), Respiratory Syncytial Virus A (RSV A)

Matrix – Transport Medium

Stability – Single use. Once thawed, use immediately

Shelf life – Up to 2 years from date of manufacture

Product Code	Product Description	Pack Size
RTXMQC01-D	Respiratory Triplex Q Control (RUO)	5 x 0.5ml (Positive) 5 x 0.5ml (Negative)

Related Products

The Qnostics range comprises Q Controls, Analytical Q Panels and Evaluation Panels for a wide range of respiratory pathogens.

- **Q Controls** - Adenovirus (ADV) | Influenza A Virus (INFA) | Influenza B Virus (INFB) | Parainfluenza Virus (PINF)
Respiratory Syncytial Virus A (RSV A) | Rhinovirus (RV)
- **Analytical Q Panels** - Influenza A Virus (INFA) | Influenza B Virus (INFB) | Parainfluenza Virus (PINF)
Respiratory Syncytial Virus A (RSV A)
- **Evaluation Panels** - Respiratory Infection Testing



Randox Laboratories Ltd

55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom | +44 (0) 28 9442 2413 | marketing@randox.com | randox.com

Information correct at time of print. Randox Laboratories Ltd is a subsidiary of Randox Holdings Limited a company registered within Northern Ireland with company number N.I. 614690. VAT Registered Number: GB 151 6827 08. Product availability may vary from country to country. Please contact your local Randox representative for information. Products may be for Research Use Only and not for use in diagnostic procedures in the USA.