





## Livatek



High flow distal Codan valve (up to 280 mL/min) and anti-siphon valves (to prevent free-flow and backflow) for the syringe pumps included



Microbore tubing used in all Tivatek™ kits ensures a low priming volume and reduces drug wastage



Anti-kink bridge as additional safety feature



High specification non-kink tubing offers a reduced risk of occlusion – even with several tight knots in the infusion line, drugs continue to run at flow rates set



2-part TIVA Infusion Sets with integral needle-free access ports and back check valve now available

## Ordering information

Total Intravenous Anaesthesia Sets - clinical dead space 0.68 mL				
Product code	Product description	Priming volume (mL)	Units per box	NHS supply chain code
2-way				
TT202E	2-way 2.0m TIVA Infusion Set with high flow valve	2.2	50	FFF903
TT302E	2-way 3.0m TIVA Infusion Set with high flow valve	3.0	50	FFF904
3-way				
TT223E	3-way 2.2m TIVA Infusion Set with high flow valve	4.3	50	FFF905
TT223ER	3-way 2.2m TIVA Infusion Set with high flow valve and red connectors	4.3	50	FFF906
TT253E	3-way 2.5m TIVA Infusion Set with high flow valve	4.9	50	FFF907
TT253FE	3-way 2.5m 2-part TIVA Infusion Set with 2 needle-free access ports and valve	5.7	50	FFF1783
TT303E	3-way 3.0m TIVA Infusion Set with high flow valve	5.7	50	FFF908
4-way				
TT224ER	4-way 2.2m TIVA Infusion Set with high flow valve and red connectors	6.4	25	FFF909
TT304E	4-way 3.0m TIVA Infusion Set with high flow valve	8.3	25	FFF910

## Manufactured in the UK











