

# Protex™ N95 Particulate Respirator Model SR9520

NIOSH Approved | TC-84A-9295

## Featuring the Protex™ All-Day-Comfort™ System

Shawmut – a world leader in advanced materials – has put everything we know into our Protex™ N95 molded-cup-style respirator, including our patent-pending Protex ADC™ All-Day-Comfort System.

Designed, engineered and manufactured in the USA, Shawmut's new Protex™ Model SR9520 mask is easier to breathe through, as the Protex™ ADC™ System addresses users' concerns that most molded-cup N95s are not comfortable to wear for long periods.

The Protex™ ADC™ System combines a unique soft but strong inner layer, incredibly lightweight but highly efficient and effective inner filtration layer, and high-sealing viscoelastic nose foam for a secure but comfortable seal.



The Protex™  
All-Day-Comfort™  
System

### Get Protex™ ADC™ for Advanced Protection and Comfort.

- A uniquely soft but strong inner layer to reduce abrasiveness on the skin
- An incredibly lightweight but highly efficient and effective inner filtration layer is easier to breathe through
- High sealing, slow-recovery viscoelastic nose foam for a secure but comfortable seal that provides better fit without excessive pressure on the face and skin
- Integrated adjustable aluminum nose clip provides a secure seal and reduced eyeglass fogging
- Expanded force curve elastic head straps are engineered to fit a wider range of sizes

**Reorder Number: SR9520 – 12**



## PRODUCT DETAILS

### PROTEX™ N95 PARTICULATE RESPIRATOR MODEL SR9520

**COLOR:** White with yellow straps

**STYLE:** Molded Cup, non-valve

**SIZE:** Medium/Large

**WEIGHT:** 0.3 oz per mask

**PACKAGING:** 20 masks/box, 12 boxes/case

**COMPLIANCE:** NIOSH 42CFR 84, TC-84A-9295

**UPC:** Box 850026398014, Case 00850026398021

**COUNTRY OF ORIGIN:** Made In USA

#### TYPE:

- Disposable Filtering Facepiece Respirator
- Filtration Efficiency
- N95 - 95% against solid particulates & non-petroleum based liquid aerosols at 0.3 microns

#### TECHNICAL DATA:

- Box Dimensions - 5 13/16" (L) X 5" (W) X 8" (H)
- Box Weight - 0.45 lb
- Case Dimensions - 18.125" L X 12.75" W X 16.25" H
- Case Weight - 6 lbs.

To learn more about the Protex™ N95 Filtering Facepiece Respirators and to place an order, please contact us at [medicalsales@shawmutcorporation.com](mailto:medicalsales@shawmutcorporation.com)



#### MATERIALS:

- Polypropylene filter and outer cover
- PET inner shell
- Elastic straps
- Aluminum nose clip
- Viscoelastic polyurethane nose foam
- All materials are latex and PVC free

#### STORAGE CONDITIONS AND SHELF LIFE:

Protex™ N95 particulate respirators have a standard shelf life of 5 years when respirators are stored in their original packaging within climatic conditions ranging from -4 °F (-20°C) to +86 °F (+30°C).

Use respirators before the “use by” date specified on each box.

#### TECHNICAL ASSISTANCE:

508-588-3300 | [shawmutcorp.com/safety](http://shawmutcorp.com/safety)