

Protex™ N95 Particulate Respirator

Model SR9520 – Size M/L

Model SR9520S – Size Small

NIOSH APPROVED | TC-84A-9295

Featuring the Protex™ All-Day-Comfort™ System



FILTER EFFICIENCIES

TYPE		EFFICIENCY	
N	Not Resistant	95	95% minimum filter efficiency in removing particulates.
R	Resistant to Oil	99	99% minimum filter efficiency in removing particulates.
P	Oil Proof	100	99.7% minimum filter efficiency in removing particulates.

Get Protex™ ADC™ for Advanced Protection and Comfort

Shawmut – a world leader in advanced materials – has put everything we know into the Protex™ N95 molded-cup-style respirator, including the patent-pending Protex™ All-Day-Comfort™ System. Designed, engineered and manufactured in the USA, Shawmut's Protex™ respirators are easier to breathe and speak through while being comfortable to wear for long periods.



The Protex™ All-Day-Comfort™ System Gives You:

- **Less breathing resistance** – proprietary filtering technology exceeds N95 standards while allowing for lighter weight and higher breathability
- **Less pressure on the face** for all-day-comfort and wear, Protex™ N95 users feel the softness of the foam with less eyeglass fogging when shaped around the bridge of the nose properly
- **Reduced moisture buildup** inside the respirator



The Protex™
All-Day-Comfort™
System

N95 Respirator Comparison Chart

Category	Protex™	Competitor A	Competitor B	Competitor C	Competitor D
100% US-Sourced Material	✓				
Country of Origin	USA	USA w/Global Materials	China	Global	India
Easy Breathing ($< 10\text{mm H}_2\text{O}$)	✓	✓		✓	
Cooler to Wear	✓				
All-Day-Comfort™ System	✓				
Fits and Seals Widest Range of Face Sizes	✓	✓			
Meets ASTM F3407 Standard	✓	✓			
Fit Capability Score	97%	73%	0%	23%	16%

Comparison results obtained through in-house laboratory testing

**Contact Shawmut Today to Order
Your N95 Respirator Made in the USA**

To Learn More or Place an Order:

MedicalSales@ShawmutCorporation.com

shawmutcorp.com/safety

Reorder Information

Part #	Size	QTY/Box	Boxes/Case
SR9520-12	M/L	20	12
SR9520S-12	Small	20	12