

# Protex™ N95 Particulate Respirator Model SR9520S (size Small)

NIOSH Approved | TC-84A-9295

## Featuring the Protex™ All-Day-Comfort™ System

Shawmut – a world leader in advanced materials – has put everything we know into our Protex™ N95 molded-cup-style respirator, including our patent-pending Protex™ All-Day-Comfort™ System.

Designed, engineered and manufactured in the USA, Shawmut's new Protex™ Model SR9520S respirator is easier to breathe through, as the Protex™ ADC™ System addresses users' concerns that most molded-cup N95s are not comfortable to wear for long periods. The small size N95 respirator provides a tight fit and secure seal on smaller faces.

The Protex™ All-Day-Comfort™ System combines a unique soft but strong inner layer, incredibly lightweight but highly efficient and effective inner filtration layer, and highsealing viscoelastic nose foam for a secure but comfortable seal.



The Protex™  
**All-Day-Comfort™**  
System

### Get Protex™ ADC™ for Advanced Protection and Comfort

Using 3D scanning, in-house testing, and rapid prototyping, Shawmut innovated the N95 Mask with its patent-pending Protex™ ADC™ all-day-comfort system to include:

- A uniquely soft but strong inner layer to reduce abrasiveness on the skin
- An incredibly lightweight but highly efficient and effective inner filtration layer is easier to breathe through
- High sealing, slow-recovery viscoelastic nose foam for a secure but comfortable seal that provides better fit without excessive pressure on the face and skin
- Integrated adjustable aluminum nose clip provides a secure seal and reduced eyeglass fogging
- Expanded force curve elastic head straps are engineered to fit a wider range of sizes



## PRODUCT DETAILS

### Protex™ N95 Particulate Respirator Model SR9520S

**Style:** Molded cup, non-valve

**Size:** Small

**Color:** White with yellow straps

**Weight:** 0.24 oz per respirator

**Packaging:** 20 N95s/box, 12 boxes/case

**Compliance:** NIOSH 42 CFR 84, TC-84A-9295

**UPC:** Box 850026398441, Case 10850026398448

**Country of Origin:** Made in USA

#### Type:

- Disposable filtering facepiece respirator
- N Series – Used in environments free of oil aerosols
- N95 – 95% against solid particulates and nonpetroleum based liquid aerosols at 0.3 microns

#### Technical Data:

- Product dimensions: 5" (L) x 4.125" (W) x 1.75" (H)
- Box dimensions: 5.8125" (L) x 5" (W) x 8" (H)
- Box weight: 0.31 lbs.
- Case Dimensions: 18.125" (L) x 12.75" (W) x 16.25" (H)
- Case Weight: 4.32 lbs.

### Contact Shawmut Today to Order Your N95 Respirators

#### To Learn More or Place an Order:

[MedicalSales@ShawmutCorporation.com](mailto:MedicalSales@ShawmutCorporation.com)  
[shawmutcorp.com/safety](http://shawmutcorp.com/safety)

**Shawmut Corporation**  
208 Manley Street  
West Bridgewater, MA 02379

© 2022 Shawmut Corporation. All Rights Reserved. INT-Lit# 00177/0 1/22



#### Materials:

- Polypropylene filter and outer cover
- PET inner shell
- Elastic straps
- Aluminum nose clip
- Viscoelastic polyurethane nose foam
- All materials are PVC free
- Not made with natural rubber latex

#### Storage Conditions and Shelf Life:

Protex™ N95 particulate respirators have a standard shelf life of 5 years when respirators are stored in their original packaging within climatic conditions ranging from -4 °F (-20 °C) to +86 °F (+30 °C).

Use respirators before the “use by” date specified on each box.

**Technical Assistance:** 508-588-3300

Reorder Information			
Part #	Size	QTY/Box	Boxes/Case
SR9520S-12	Small	20	12