# Protex™ N95 Particulate Respirator Model SR9520 (size M/L)

NIOSH Approved TC-84A-9295

### Featuring the Protex™ All-Day-Comfort™ System

Made by Shawmut Corporation, a world leader in advanced materials, the Protex™ SR9520 is a NIOSH approved N95 disposable respirator and meets NIOSH 42 CFR 84 N95 requirements for a minimum 95% filtration efficiency.

Designed, engineered, and manufactured in America, the cup-style Protex™ SR9520 respirator was created with patent pending ADC Protex All-Day-Comfort System® technology so you can count on it to fit better, protect better, and breathe easier all day long.

An N95 respirator that doesn't fit properly may not protect you. Protex™ N95s are designed to fit the widest range of face sizes and shapes and have received Fit Capability scores of up to 97% when tested against the ASTM F3407-20 Respirator Fit Capability (RFC) standard.





## Get Protex<sup>™</sup> ADC<sup>™</sup> for Advanced Protection and Comfort

The Protex<sup>™</sup> ADC<sup>™</sup> combines a unique soft but strong inner layer, incredibly lightweight but highly efficient and effective inner filtration layer, high-sealing viscoelastic nose foam, and expanded force curve elastic head straps for a secure but comfortable seal.

- MORE SECURE FIT Protex<sup>™</sup> N95s are designed to fit a wider range of face shapes and sizes. And a better fit means better performance: Protex<sup>™</sup> N95s remove more than 95% of all 0.3+ micron particles.
- BEST-IN-CLASS BREATHABILITY –
   Proprietary Shawmut MB™ filtering technology exceeds NIOSH N95 standards while delivering better breathability and cooler, fresher air inside the cup.
- ALL-DAY COMFORT Protex<sup>™</sup> N95s feature softer, flanged edges that provide a more comfortable fit with less pressure on your face.
- LESS EYEGLASS FOGGING When formed around your nose and worn properly, the unique viscoelastic nose foam and aluminum nose clip reduce eyeglass fogging for better visibility.





# shawmut

### **PRODUCT DETAILS**

Protex<sup>™</sup> N95 Particulate Respirator Model SR9520

Style: Molded cup, non-valve

Size: Medium/Large

**Color:** White with yellow straps

Weight: 0.3 oz per respirator

Packaging: 20 N95s/box, 12 boxes/case

**Compliance:** NIOSH 42 CFR 84, TC-84A-9295

**UPC:** Box 850026398014, Case 00850026398021

Country of Origin: Made in USA

#### Type:

- Disposable filtering facepiece respirator
- N Series Used in environments free of oil aerosols
- N95 95% against solid particulates and nonpetroleum based liquid aerosols at 0.3 microns

#### **Technical Data:**

- Box dimensions: 5.8125" (L) x 5" (W) x 8" (H)
- Box weight: 0.45 lbs.
- Case Dimensions: 18.125" (L) x 12.75" (W) x 16.25" (H)
- · Case Weight: 6 lbs.

### Contact Shawmut Today to Order Your Protex™ N95 Respirators

#### To Learn More or Place an Order:

<u>PPESales@ShawmutCorporation.com</u> shawmutcorp.com/safety

**Shawmut Corporation** 208 Manley Street West Bridgewater, MA 02379



#### Materials:

- Polypropylene filter and outer cover
- PET inner shell
- Elastic straps
- Aluminum nose clip
- · Viscoelastic polyurethane nose foam
- · All materials are PVC free
- Not made with natural rubber latex

#### **Storage Conditions and Shelf Life:**

Protex<sup>™</sup> N95 particulate respirators have a standard shelf life of 5 years when respirators are stored in their original packaging within climatic conditions ranging from -4  $^{\circ}$ F (-20  $^{\circ}$ C) to +86  $^{\circ}$ F (+30  $^{\circ}$ C).

Use respirators before the "use by" date specified on each box.

Technical Assistance: 508-588-3300

Reorder Information			
Part #	Size	QTY/Box	Boxes/Case
SR9520-12	M/L	20	12

