

Esophageal Stents

Gastroenterology



ALIMAXX-ES™
Fully Covered Esophageal Stent

EndoMAXX®
Fully Covered Esophageal Stent



Esophageal Stents

CORE ESOPHAGEAL STENT TECHNOLOGY


Anti-Migration Struts

Designed to reduce stent migration

Laser Cut Design

- Virtually non-foreshortening
- Improved deployment accuracy

Advanced Patient-Focused Design



SOFTER

FIRMER

SOFTER

Single User Deployment System

- One-handed deployment mechanism
- Second hand free to control stent position

EndoMAXX®

Fully Covered Esophageal Stent

Unique Features

- Migration rates as low as 11%*
- Distal and Proximal metal sutures




ALIMAXX-ES™

Fully Covered Esophageal Stent

Unique Features

- Expansion force as low as 1.6 Newtons, compared to 5.4 Newtons found in woven stents**
- Small diameter sizes



* Dua et al. GIE. Efficacy and safety of a new fully covered self-expandable non-foreshortening metal esophageal stent. 2014.

** Data on file

Esophageal Stents

ENDOMAXX® FULLY COVERED ESOPHAGEAL STENT



Catalog Number	Stent Size Options			Delivery Device		Packaged
	Stent Body Diameter (mm)	Expanded Stent Length (mm)	Stent Flare Diameter (mm)	Working Diameter (F)	Working Length (cm)	

OVER-THE-WIRE

MAXX-197	19	70	24	24	90	1 unit
MAXX-1910	19	100	24	24	90	1 unit
MAXX-1912	19	120	24	24	90	1 unit
MAXX-1915	19	150	24	24	90	1 unit
MAXX-237	23	70	28	24	90	1 unit
MAXX-2310	23	100	28	24	90	1 unit
MAXX-2312	23	120	28	24	90	1 unit
MAXX-2315	23	150	28	24	90	1 unit

UNIQUE SIZES are noted with **RED/BOLD** part numbers

ALIMAXX-ES™ FULLY COVERED ESOPHAGEAL STENT



Catalog Number	Stent Size Options		Delivery Device		Packaged
	Nominal Stent Diameter (mm)	Expanded Stent Diameter (mm)	Working Diameter (F)	Working Length (cm)	

OVER-THE-WIRE

80129-207	12	70	16	67	1 unit
80129-208	12	100	16	64	1 unit
80129-209	12	120	16	62	1 unit
80129-210	14	70	16	67	1 unit
80129-211	14	100	16	64	1 unit
80129-212	14	120	16	62	1 unit
80129-213	16	70	16	67	1 unit
80129-214	16	100	16	64	1 unit
80129-215	16	120	16	62	1 unit
80129-201	18	70	16	67	1 unit
80129-202	18	100	16	64	1 unit
80129-203	18	120	16	62	1 unit
80129-204	22	70	16	67	1 unit
80129-205	22	100	16	64	1 unit
80129-206	22	120	16	62	1 unit

ADDITIONAL STRICTURE MANAGEMENT PRODUCTS



BIG60®
60 mL inflation device

Elation5™
5-stage dilation balloon



Elation®
3-stage dilation balloon



Elation®
Wire guided dilation balloon



Maxxwire®
Guidewire designed for stent placement

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.