



CRYOGENIC ACCESSORIES

CRYO-APRON

Cryogenic apron is waterproof and offers superior thermal protection. Select the appropriate length for your specific application.

Highlighted Features

- Superior thermal protection
- Waterproof
- Comfortable / Lightweight / Breathable
- Warmth and splash protection for torso and legs
- Adjustable straps at neck and waist
- Side release buckles for quick on/off

Typical Applications

- Cold rooms / walk-in freezers
- Blast freezers
- Delivery and filling of cryogenic liquids
- Hazardous material cryo-liquid spill response

Environments

- Cryo preservation facilities
- Cryogenic liquid producers
- Clinical laboratory
- Bio-Medical
- Blood Banks
- Frozen food processing
- Pharmaceutical



SAFETY FIRST
WEAR PROTECTIVE APPAREL
WHEN RISK IS PRESENT

**KEEP GLOVES AND APRONS
READILY ACCESSIBLE**

*Caution: Not Intended for Immersion in
Liquid Nitrogen or other Cryogenic Liquids*



Selection & Usage

- Proper fit contributes to product performance – a tight fit leads to thermal loss.
- Select the appropriate glove length and apron length for your specific application.
- Select **Waterproof Cryo-Gloves** for use in cryogenic atmospheres – glove fit should be loose enough to allow for quick removal if necessary.
- Select **Cryo-Industrial Gloves** when gloves are in contact with abrasive surfaces.
- Store gloves and aprons in a clean dry space.
- Periodically inspect the condition of your gloves and aprons – replace garments that show excessive wear.

Order Information

Order Catalog Number	Description	Length (in.)	Length (cm)	Estimated Weight
CRY AP36	Cryo-Apron 36" long	36	90.85	0.64 lbs.
CRY AP42	Cryo-Apron 42" long	42	100.4	0.64 lbs.
CRY AP48	Cryo-Apron 48" long	48	121.75	0.64 lbs.
CRY AP54	Cryo-Apron 54" long	54	137.0	0.64 lbs.

Can't find what you need? We will be happy to assist you to find the perfect personal protective equipment to fit your needs.