



CRYOGENIC GLOVES

14" to 15" Mid-Arm Length Standard Cryogenic Gloves

Mid-Arm length, approximately 14" to 15" standard cryogenic gloves.

Highlighted Features

- Superior thermal protection
- Comfortable / Lightweight / Excellent Dexterity
- Breathable / Stay-Dry lining
- Protection from splashes of cryogenic liquids
- Extra comfort when working in cryogenic atmospheres for extended periods

Typical Applications

- LN2 atmospheres and environments
- Refrigerated liquid baths
- Cryo freezers
- Cryogenic grinding

Environments

- Fertility Clinics
- Bio-Medical
- Industrial
- Pharmaceutical
- Genetic and Bio Technology



SAFETY FIRST
WEAR PROTECTIVE APPAREL
WHEN RISK IS PRESENT

KEEP GLOVES AND APRONS
READILY ACCESSIBLE

*Caution: Not Intended for Immersion in
Liquid Nitrogen or other Cryogenic Liquids*

Selection & Usage

- Proper fit contributes to product performance – a tight fit leads to thermal loss.
- Select the appropriate glove length and apron length for your specific application.
- Select **Waterproof Cryo-Gloves** for use in cryogenic atmospheres – glove fit should be loose enough to allow for quick removal if necessary.
- Select **Cryo-Industrial Gloves** when gloves are in contact with abrasive surfaces.
- Store gloves and aprons in a clean dry space.
- Periodically inspect the condition of your gloves and aprons – replace garments that show excessive wear.



Order Information

Mid-Arm Length – 13 1/2" - 15 1/4" (345-390 mm)

Order Catalog Number	Description	Size	Estimated Weight
CRY-GL1415-XL	Mid-Arm length cryogenic glove (1 pair)	Extra Large	0.62 lbs
CRY-GL1415-L	Mid-Arm length cryogenic glove (1 pair)	Large	0.62 lbs
CRY-GL1415-M	Mid-Arm length cryogenic glove (1 pair)	Medium	0.62 lbs
CRY-GL1415-S	Mid-Arm length cryogenic glove (1 pair)	Small	0.62 lbs

Can't find what you need? We will be happy to assist you to find the perfect cryogenic glove to fit your needs.