



# CRYOGENIC GLOVES

25" to 27" Shoulder Length  
Standard Cryogenic Gloves

Shoulder length, approximately 25" to 27" standard cryogenic gloves.

## Highlighted Features

- Superior thermal protection
- Comfortable / Lightweight / Excellent Dexterity
- Breathable / Stay-Dry lining
- Protection from splashes of cryogenic liquids
- Extra comfort when working in cryogenic atmospheres for extended periods

## Typical Applications

- LN2 atmospheres and environments
- Refrigerated liquid baths
- Cryo freezers
- Cryogenic grinding

## Environments

- Fertility Clinics
- Bio-Medical
- Industrial
- Pharmaceutical
- Genetic and Bio Technology



**SAFETY FIRST**  
WEAR PROTECTIVE APPAREL  
WHEN RISK IS PRESENT

**KEEP GLOVES AND APRONS  
READILY ACCESSIBLE**

*Caution: Not Intended for Immersion in  
Liquid Nitrogen or other Cryogenic Liquids*

### Selection & Usage

- Proper fit contributes to product performance – a tight fit leads to thermal loss.
- Select the appropriate glove length and apron length for your specific application.
- Select **Waterproof Cryo-Gloves** for use in cryogenic atmospheres – glove fit should be loose enough to allow for quick removal if necessary.
- Select **Cryo-Industrial Gloves** when gloves are in contact with abrasive surfaces.
- Store gloves and aprons in a clean dry space.
- Periodically inspect the condition of your gloves and aprons – replace garments that show excessive wear.



## Order Information

Shoulder Length – 24 3/4" - 27" (625-685 mm)

Order Catalog Number	Description	Size	Estimated Weight
CRY-GL2527-XL	Shoulder length cryogenic glove (1 pair)	Extra Large	0.62 lbs
CRY- GL2527-L	Shoulder length cryogenic glove (1 pair)	Large	0.62 lbs
CRY- GL2527-M	Shoulder length cryogenic glove (1 pair)	Medium	0.62 lbs
CRY- GL2527-S	Shoulder length cryogenic glove (1 pair)	Small	0.62 lbs

*Can't find what you need? We will be happy to assist you to find the perfect cryogenic glove to fit your needs.*