




MOSAIC
SAMPLE MANAGEMENT SOFTWARE

Mosaic 8 Integrations



Mosaic software helps to digitalise your lab by offering drivers and support for tight 'driven' integrations with an extensive range of instruments and devices from multiple vendors. Auto-validated file integration is possible with any lab automation



MOSAIC SOFTWARE SUPPORTS THE FOLLOWING

Automated Stores

- Brooks Matrical MiniStore
- Brooks REMP LSS/MSS
- Brooks REMP SSS
- Brooks RTS Flexi-Store
- Brooks SampleStore™ II with Automation Integration Module (AIM)
- Brooks TAP Haystack
- Hamilton ASM
- Hamilton BiOS
- Hamilton SAM HD
- Hamilton Verso
- Modula Store
- NEXUS Biosystems Universal Store
- SPT Labtech comPOUND®
- SPT Labtech arctic®

Manual Stores

- Kardex Shuttle

Biological Registration

- Dassault Systemes BIOVIA™ Registration
- Dotmatics BioRegister
- Genedata Biologics
- Watson LIMS

Screening Data Analysis

- Dotmatics Ultra Screening
- Genedata Screener
- IDBS ActivityBase
- PerkinElmer Signals Screening

Small Molecule Registration

- CDD Vault Registration System
- ChemAxon Chemical Registration
- Dassault Systemes BIOVIA™ Chemical Registration
- Deltasoft ChemCart
- Dotmatics Register
- IDBS ActivityBase
- Scilligence RegMol

Other Software

- Certara D360
- ChemAxon Compliance Checker
- Dassault Systemes BIOVIA™ Workbook
- Dotmatics Browser

Dry Powder Handling

- Mettler-Toledo balances
- Mettler-Toledo FlexiWeigh®
- Sartorius balances

Liquid Handling

- Agilent BioCel
- Agilent BenchCel
- Agilent Bravo®
- Agilent VPrep®
- Beckman Coulter Access™ Dual Robot System
- Beckman Coulter Access™ Laboratory Workstation
- Beckman Coulter Biomek® FX
- Beckman Coulter Biomek® NX
- Beckman Coulter Echo® standalone (formerly Labcyte)



- Beckman Coulter POD™
- Hamilton STAR
- Hamilton STARlet
- Hamilton Microlab STAR
- Hamilton Microlab STARplus
- Hamilton Microlab VANTAGE
- HighRes Biosolutions ACell
- HighRes Biosolutions CoLab Flex
- HighRes Biosolutions MicroStar
- HighRes Biosolutions NanoCell
- MatriCal AutoMAP
- Perkin Elmer MultiPROBE II
- SPT Labtech mosquito®
- Tecan Freedom EVO®
- Tecan Fluent®
- Zinsser LISSY®
- Biosero BioRead A6 High Speed Single Rack Scanner
- Biosero BioRead AV6 Deep Focus Single Rack Scanner
- Biosero BioRead Express Rack Reader
- Biosero BioRead Multirack Scanner
- Biosero DataPac Cube Scanner
- FluidX Impression™ Rapid Rack Scanner
- FluidX Perception™ Rack Reader
- FluidX Xtr-96 2D Barcode Reader
- Micronic Tracxer Code Reader
- SPT Labtech cube Rack Scanner
- Thermo Scientific VisionMate™ 96 2D Barcode Reader
- Thermo Scientific VisionMate™ HS 2D Barcode Reader
- Thermo Scientific VisionMate™ SR 2D Barcode Reader
- Ziath DataPac™ Cube Rack Reader
- Ziath DataPac™ Deep Focus Rack Reader
- Ziath DataPac™ Express Rack Reader
- Ziath DataPac™ High Speed Reader
- Ziath DataPac™ Multirack Reader
- Ziath DataPac™ Single Tube Reader

2D Tube Sorting

- BioMicroLab XL9
- BioMicroLab XL20
- BioMicroLab XL100
- BioMicroLab XL200
- Ziath Mohawk

2D Tube Scanning

- AltemisLab SampleScan Rack Barcode Reader
- BioMicroLab SampleScan™ Mini 2D Barcode Reader
- BioMicroLab SampleScan™ Mini Automated 2D Barcode Reader
- BioMicroLab SampleScan™ 96 Barcode Reader
- BioMicroLab SampleScan™ 96 Automated Barcode Reader
- BioMicroLab SampleScan™ HS Barcode Reader
- BioMicroLab SampleScan™ E-Series Barcode Reader
- BioMicroLab SampleScan™ E-Series Automated Barcode Reader

Other

- Any “keyboard wedge” 2D barcode reader
- Any label printer via BarTender labelling
- Agilent VCode labeller
- Agilent BenchCel labeller
- BioMicroLab VolumeCheck™
- Brooks Tube Auditor™
- BarTender labelling
- CyBio® CyBi-Print labeller
- LogistiVIEW Digital Workforce Platform
- Vuzix Smart Glasses M300/M300XL/M400





Different device?

Almost any lab equipment can be integrated with Mosaic software using our validated file exchange or JFM Application Programming Interfaces.

If you don't see the device you want in this document, please talk to a Titian representative.

Thank you for your interest in this document



White papers



Application notes



Product brochures



Video content



Webinar recordings

You can find a collection of further related material in the resources section on our website.

info@titian.co.uk

www.titian.co.uk

UK Tel:

+44 20 7367 6869

USA Tel:

+1 508 366 2234