EMBODIED LABS™ SETUP MODES

Modes	What is it?	Why use this mode?
Individual in the VR Headset Setup Guide	One person wears the VR headset at a time. They sit on a swivel chair to go through the full 360 degrees experience.	 This mode is the most immersive: Learners feel emotionally connected with they see, hear, and feel. Retention. Learners remember this experience for a long time.
VR with an In-Person Group Setup Guide	One person wears the VR headset while a group in the same room follows along on a TV screen connected to the headset.	This mode offers the opportunity for a group of learners to debrief together after each lab section, to share observations and experiences, even personal stories, to deepen their learning.
VR with a group over Zoom Setup Guide	An instructor shares their in-VR headset experience with a group over Zoom	This shares the Instructor's live, fully immersive VR experience over Zoom, so the instructor can choose and control how the experience unfolds. The same debrief discussion opportunity exists as with VR to an in-person group.
Embodied Labs Online Setup Guide	 Log in to app.embodiedlabs.com on Google Chrome browser Learners navigate through the lab experiences by clicking and dragging with their mouse. No headset is needed. 	 This mode offers flexibility: Use it when there isn't time for each learner do an individual in-headset experience so they do part of a lab in the headset and the rest in online mode. Use it to give learners a chance to further review on their own time a lab experience that they already did once as a group.
Embodied Labs Online Remote Instruction Setup Guide	An instructor shares their screen on Microsoft Teams or Zoom to lead a group of learners through a lab in online format.	Ease of use. The instructor needs little equipment and can fully choose and control how the experience unfolds based on desired learning outcomes. Direct follow up is easy too, by at the end of a session assigning learners to log into app.embodiedlabs.com and do additional modules or labs or to review the same lab again on their own time.

